

APPLICATION NOTE: SUSTAINABLE SOIL MANAGEMENT

Techniques

- NIR Spectroscopy
- Field Spectroscopy

Keywords

- Smart Agriculture
- Bulk Soil Density
- NIR Diffuse Reflection Soil Moisture Content
- Soil Classification
- Sustainable farming

Introduction

Agricultural research of today goes into feeding the world's population tomorrow. Year after year, farmers try to produce more food out of the same resources, and rapidly changing climate conditions only heighten the pressure on the world's food supply. Changing weather patterns means droughts deepen, floods are more frequent, extremes of heat and cold swing wider, and so inevitably we must produce more with less. Shorter growing seasons, fewer resources such as clean water, and even the soil itself can become depleted. Right



Layers of soil showing compacted zone

now, scientists are hard at work developing the tools and technologies that will make the future of farming possible. Spectroscopy is a key enabling technology for many of those tools and technologies in development today. From innovative research to integration in sensors and analytics, spectroscopy is everywhere. And so is Avantes. Trusted in labs, production lines and field research outposts, Avantes instruments deliver proven results around the world.

Soil Management

Soil is a compound mixture of organic matter, minerals, gases, liquids, and even living organisms. In addition to supporting the plant life we need for crops, soil also functions as a means of storing, transpor-



Irrigation system

ting, and purifying water; it helps to modify the atmosphere we all depend on, and even serves as habitat for organisms large and small. Sustainable soil management is as critical for future food production as it is for life on Earth.

Soil health is elemental to sustainable land management and is an important consideration for farms of all sizes. Anything from erosion to contamination, loss of biodiversity, soil compaction, and everything in between, can be detrimental to crop production and the viability of the farm itself.

Numerous studies and technologies are in development for analyzing and managing soil health and Avantes is at the forefront of research and technology development, protecting our future food supply.

Moisture Measurements

Historically, measurements of soil moisture have employed a device called a tensiometer which uses a hollow tube with a porous reservoir of water on top and a gauge. The tube is inserted into the soil and the water is drawn into the soil from the cup, creating a vacuum, until it reaches equilibrium. The gauge allows the user to gather a reading based on that vacuum that correlates to water carrying potential in the soil matrix. This data allows farmers

> Avantes BV Oude Apeldoornseweg 28 7333 NS Apeldoorn The Netherlands Tel. +31 313 670 170 info@avantes.com www.avantes.com

Avantes USA 500 S. Arthur Ave., Unit 500 Louisville. CO 80027 Tel. +1 (303)-410-8668 infousa@avantes.com www.avantesusa.com



MEMBER OF THE NYNOMIC GROUP

to determine the need for irrigation, but there are drawbacks. The tensiometer is slow, requiring a period of time for water to reach equilibrium which limits the scale of use.

As early as 1978, researchers have been interested in the effects of soil moisture on, first the visible spectrum, but as spectroscopy techniques improves, eventually extending into the near- and mid-infrared. [2] In 2014, researchers in Hungary were working to calibrate spectral data intended to develop algorithms that would allow fast, field-scale measurements of soil moisture content [3].

The development of algorithms that will one day allow for field-scale deployment of spectroscopy-based moisture measurements begins first in the laboratory with data collection. Soil samples collected from orchards with varying soil characteristics from around the region were first kiln dried in the lab to a consistent aridity. Water was then reintroduced slowly, with spectra collected each 2.5ml until the samples were fully saturated. They identified the wavelengths 1450-1460 nm and 1920-1930 nm as the most sensitive for quantifying soil moisture. The AvaSpec-NIR256-2.5-HSC-EVO and AvaSpec-NIR256-1.7-EVO allow the range and sensitivity for this type of application in the laboratory, but future developments in compact

Effects of compacted Soil

NIR spectrometers, such as the AvaSpec-CompactLine NIR, expected in 2020, could one day lead to replacing the tensiometer with handheld field instruments for rapid assessment of moisture content and integration with irrigation systems and to calibrate and authenticate (ground truth) airborne hyperspectral imaging technology. [3]

Soil Characterization

Soil is not a homogenous mixture, differences in minerals, organic matter,

and particle size, as just a few examples, can alter any number of soil characteristics. [4] Soil type will determine everything from irrigation schedules to types of crops likely to prosper. The US Department of Agriculture recognizes 12 types of soil with sand at one end of the spectrum and clay at the other. [5] Each type has predictable characteristics, including color. In fact, the most common

method for determining soil type relies on subjective (and fallible) personal experience to compare against a specially designed color chart from the Munsell Corporation. [2]

Other methods for soil classification require chemical processes that might have adverse effects on data interpretation. [3] They are also time-consuming and require more advanced technical skills to perform. Optical spectral sampling, on the other

> hand, requires little or no sample preparation and no harsh chemicals, and several parameters can be analyzed from the same spectral data.

Researchers in Hamadan, Iran investigated UV/VIS/NIR spectroscopy to analyze a number of soil parameters including color, pH, electrical conductivity, moisture content, available organic carbon, total nitrogen, and exchangeable cations, as well

as identifying minerals such as iron, titanium oxides, calcium, magnesium, potassium, and sodium, just to name a few. [3] Samples were collected, randomized, and kiln-dried then ground and sieved before being divided into calibration and validation sets. A full battery of chemometric data was collected from the set of validation samples.

An average of more than 24 broadband spectral scans were collected per sample in the calibration set using the AvaSpec-



Soil samples for characterization testing

ULS2048 spectrometer to cover the UV-VIS range from 200-1100 nm at 1 nm resolution. For the 1000-2400 nm range, the AvaSpec- AvaSpec-NIR256-2.5-HSC-EVO allows measurements in the NIR. Spectral analysis performed on the calibration set allowed a series of operations for statistical analysis based on principal component analysis (PCA) and least partial squares regression analysis (LPS-R) to relate spectral reflectance measurements to observable soil properties. Based on their results, these scientists determined that they could sufficiently validate soil classification using a fraction of available variables.

Bulk Density and Soil Compression

Soil compaction is the compression of spaces between soil particles that would normally hold air or water. It leads to poor root development, oxygen deficiency, and other deficiencies. These factors ultimately reduce crop quality and yields. Compaction can be man-made or the result of natural processes, but is a serious environmental and agricultural problem. [1] Common causes include extensive use of heavy machinery, repeated plowing to

system. [1]

the measurement of gravimetric moisture

content, to allow calculation of bulk den-

sity (BD) with one easy-to-manage probe

This novel probe design combined a

penetrometer for the measurement of

soil penetration resistance, with a hol-

low shaft housing optical fibers coupled

to a sapphire window in the shaft body

at one end, and on the other, to the 20

watt halogen lamp for input, and to the

AvaSpec-NIR256-2.5 spectrometer for NIR

reflection data output and capture. Time

 (Θv) achieved by integrating an electrode

in the shaft of the probe in the form of

a copper ring insulated from the probe

body. Hundreds of measurements were

(ANN) to model the fusion of data and

used to train an artificial neural network

domain reflectometry measurements

quantify volumetric moisture content



MEMBER OF THE NYNOMIC GROUP



NIR Measurements in the lab uniform depth, and use by large animal populations. Sustainable agricultural systems must work to manage soil compaction beginning with the measurement of key parameters associated with soil compaction.

Bulk soil density (BD), representing the ratio of volumetric moisture (Θ v) content to gravimetric moisture (ω) [expressed as BD = Θ v/ ω], is positively correlated with soil compaction, i.e. the higher the soil density, the greater the degree of

Avantes Meets the Challenge

The need for affordable, efficient, and reliable data is not without challenges. Many available systems are not suitable for field deployment. Data collection can also be compromised because of operator errors and inconstancies. Much of the work of designing all-in-one spectroscopybased solutions appears to focus on eliminating opportunities for user error or inconsistencies in deployment. [2] Working with an experienced partner in compaction. [1] Current technology for measuring soil compaction is the penetrometer, an instrument that, at its heart, hasn't advanced a great deal beyond pointy stick. The penetrometer is meant to mimic a growing plant root and involves a pressure gauge atop a graduated driving shaft tipped with a wider

30-degree stainless steel cone.

Research published in 2018 in the journal Computers and Electronics in Agriculture unveiled a prototype measuring system that combines a traditional penetrometer with fiber-coupled NIRS sensor. [1] The system operated in situ to combine measurement of soil penetration resistance, frequency domain reflectometry to analyze volumetric moisture content, and near-infrared diffuse reflectance for

spectroscopy can help the system design engineer to discover solutions to any number of field deployment challenges. Avantes brings 25 years of experience working hand-in-hand with researchers and equipment manufacturers to design systems that meet and exceed measurement requirements. Test drive an Avantes instrument with our exclusive demo program and discover the Avantes advantage today.



AvaSpec-NIR256/512-2.5-HSC-EVO (TEC Cooled)

References

- [1] Al-Asadi, Raed A., and Abdul M. Mouazen. "A Prototype Measuring System of Soil Bulk Density with Combined Frequency Domain Reflectometry and Visible and near Infrared Spectroscopy." Computers and Electronics in Agriculture. Elsevier, 30 June 2018. Web. 19 Nov. 2019. https://www.sciencedirect.com/science/article/pii/ S0168169917316277
- [2] Han, Pengcheng, et al. "A Smartphone-Based Soil Color Sensor: For Soil Type Classification." Computers and Electronics in Agriculture, vol. 123, 2016, pp. 232–241., doi:10.1016/j.compag.2016.02.024. https://www.sciencedirect.com/science/article/pii/ S0168169916300618
- [3] Monavar, Hosna Mohamdi. "34th International Con-

ference on Food and Agricultural Engineering." The Ires, Determination of Several Soil Properties Based on Ultra-Violet, Visible, and near-Infrared Reflectance Spectroscopy, 12 Mar. 2016, https://www.researchgate. net/publication/297322303_Determination_of_several_soil_properties_based_on_ultra-violet_visible_and_ near-infrared_reflectance_spectroscopy

- [4] "Soil." Wikipedia, Wikimedia Foundation, <u>https://en.wikipedia.org/wiki/Soil</u>.
- [5] USDA. "Natural Resources Conservation Service." The Color of Soil | NRCS Soils, US Department of Agriculture, <u>https://www.nrcs.usda.gov/wps/portal/nrcs/detail/</u> <u>soils/edu/?cid=nrcs142p2_054286</u>